

HBB-M

Technical data

Dimensions: 660 x 655 x 1,100 mm

Weight: Approx. 110 kg

Power supply and consumption:

Voltage: 400 V, 50 Hz, 3-phase, or as required

Neutral conductor: Not required

Power consumption: Approx. 2 kVA

Electrical switch cabinet:

Programmable controller: Simatic S7

Protection class: IP 44

Insulation class: B

Compressed air supply and consumption:

Pressure setting: 6 bar

Sample material which can be processed:

Sample material: Slag, sinter and ore

Body diagonal: 15 mm max.

Attainable granuable size: Approx. 1 mm to 8 mm depending on the material

Temperature: Max. 500 °C

The design of the machine complies with the applicable accident prevention and VDE regulations.

We reserve the right to make technical modifications.