

HSM 50

Technical data

Color:	white/blue
Labeling text:	English
Instruction manual:	1 copy, English
Accessories:	1 set of wrenches

Dimensions

Machine including packing:	850 x 850 x 1,100 mm
-----------------------------------	----------------------

Weight

Machine including packing:	180 kg
-----------------------------------	--------

Power supply and consumption:

Voltage:	400 V, 50 Hz, 3 Phases
Neutral conductor:	Not required
Connected load:	2,5 kVA

Electrical control cabinet (integrated):

Programmable controller:	Simatic C7
Control voltage:	24 V
Protection class:	IP 44
Insulation class:	B

Options:

- Grinding vessel 10 ccm (chrome steel or tungsten carbide)
- Grinding vessel 50 ccm (chrome steel or tungsten carbide)

The design of the machine complies with the applicable accident prevention and VDE regulations.

We reserve the right to make technical modifications.